



∞
B I O S



BIOS

BIOS collection -from the Greek word “life”- is inspired by the old apothecaries. It is made with aluminium and paper, biodegradable materials that help to preserve the environment. Cosmetics have a fresh lemon fragrance and include a paper-wrapped soap without plastic.

BIOS COSMETICS

- Aluminium tubes:

- Bath Gel (1.01 fl.oz /30ml)
- Shampoo (1.01 fl.oz /30ml)
- Conditioner (1.01 fl.oz /30ml)
- Body Milk (1.01 fl.oz /30ml)

- Paper-Wrapped Soap (1.23 oz /35gr)



BIOS

BIOS paper flowpacks are the ecological solution adapted to the most classic format. Its PSM material and its simple design make this collection fit in any establishment.

BIOS COMPLEMENTS

Paper flowpacks complements:

- Tissues (3 units).
- Biocompostable Shower Cap.
- Vanity Kit.



mcamp@armony-usa.com · T. 786.431.8387 · www.armony-usa.com



BIOS

Taking care of your guests and protecting the environment is possible with the **BIOS** dispensers: using them reduce the quantity of waste and CO2 in the atmosphere without refuse to natural cosmetic.

- **BIOS** collection has lemon fragrance, a soft formula perfect for men and women.
- Tough and long-lasting holder.
- Hermetic system: don't allow stealing or product handling.
- Gallons available to refill product.
- Bottle capacity: 9 fl.oz / 266ml
- Bottles made of 100% recyclable PET plastic.
- **BIOS** has available small bottles and paper flowpack accesories.



TYPES OF HOLDERS AVAILABLES:

- Single holder.
- Double holder.
- Triple holder.



mcamp@armony-usa.com · T. 786.431.8387 · www.armony-usa.com



BIOS

BIOS GALLONS AVAILABLE

BIOS BRAND:

- Hand & body soap gallon.
- Shampoo gallon.
- Conditioner gallon.
- Body milk gallon.



BIOS ESSENTIAL BRAND:

- Hand & body soap Essential.
- Shampoo Essential.



BIOS

TECHNICAL SHEET AND ASSEMBLY

Stainless steel bracket with security key:

- Moisture resistant.
- Cannot be tampered with.
- Screwed to the wall.
- Single holder: four anchorage points.
- Double holder: four anchorage points.
- Triple holder: four anchorage points.

Method of use:

- 1 - Anchor the holder to the wall.
- 2 - Introduce the security key into the slot and turn it ninety degrees to the left.
- 3 - Remove the metal trim.
- 4 - Place the neck of the bottles in the groove of the holder in the desired order.
- 5 - Replace the bottle holder.
- 6 - Introduce the security key in the slot and turn it ninety degrees to the right to lock the holder.

